



Column:

MEGA-DEX B-SE - 0.25mm, 0.25 μ m, 30m

Catalog Code: DEX-BSE-025-025-30

Conditions:

Injection: Split 230°C, 1:50 split ratio.

Sample Dilution: 1mg/1mL each compound.

Detector: FID 230°C.

Oven Program: 40°C, 2°C/min.

Carrier Gas: Hydrogen, 80kPa.

1. Ethyl 2-methylbutyrate (+/-)
2. 2-methylbutylacetate (+/-)
3. 2-methylbutanol (+/-)
4. 2-methylbutyrate (+/-)
5. Carvone (+/-)



since
1980



MEGA[®]
improve your GC analysis

contact us: info@mega.mi.it



follow us [@MEGAColumns](https://twitter.com/MEGAColumns)

www.mega.mi.it